

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

High Resolution Chest

Indications :

Scan Type	Helical
Scan Start\ End Location	Above apices through the tops of the kidneys
Scan Details	Patient is scanned supine and prone 1mm X 10mm slices in lung and soft tissue window
IV Contrast	N/A
Scan Delay	N/A
Reformat Details	3mm X 3mm Coronals and sagittals on supine position in lung window.
Comments	If patient is unable to lay prone than only supine images are done.
Images to send to PACS	Scouts, axials, coronals, and sagittals