

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Routine Abdomen

Indications : pain, nausea, vomiting, constipation, diarrhea,

Scan Type	Helical
Scan Start\ End Location	Above diaphragm to iliac crest
Scan Details	5mm X 5mm in abdomen window
IV Contrast	70ml injected at 2cc/sec
Scan Delay	65 second delay (venous phase)
Reformat Details	Coronals in abdomen window
Comments	Patient drinks oral contrast for this exam
Images to send to PACS	Scouts, axials and coronals