BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Stone Protocol

Indications: stone, gross hematuria, micro hematuria, flank pain, urinary frequency, hx of stones

Scan Type	Helical
Scan Start\ End Location	Diaphragm to below lesser trochanters
Scan Details	5mm X 5mm in abdomen window Non contrast
IV Contrast	N/A
Scan Delay	N/A
Reformat Details	Coronal images in abdomen window
Comments	Always done without oral and IV contrast
Images to send to PACS	Scouts, axials, and coronals in soft tissue window